


BD Veritor™ System

# Connectivity Starter Packs

Simplify Workflow.  
Enable Integration.



Deliver fast, reliable respiratory results—and ensure they're documented accurately and automatically. BD Veritor™ Connectivity Starter Packs make it easy to order tests and connectivity together by reducing manual steps.



## Two ways to get connected

### BD Veritor™ Total Connectivity Starter Pack

Tests + InfoWiFi Module + Analyzer

Best for new placements or sites expanding point of care testing.

### BD Veritor™ Connectivity Expansion Bundle

Tests + InfoWiFi Module

Designed for existing BD Veritor™ Plus Analyzer users looking to automate result reporting.



Get Started at [bdveritor.com](https://bdveritor.com)



# Why Connectivity Matters




Connectivity helps results move seamlessly from testing to documentation—so care teams can focus on patients, not paperwork.

- Reduces transcription and documentation errors
- Saves staff time and administrative effort
- Enables real time access to results


## BD Veritor™ System at a Glance



- Digitally read results in 15 minutes or less\*



- Flexible workflows: Walk Away or Analyze Now



- Clear positive/negative display



- Compact, portable design



- Single button operation, no monthly calibration

## How Connectivity Works

Operator ID



Specimen ID



Lot & Exp



BD Veritor™ Plus Analyzer equipped with InfoWifi module



Connectivity application deployed within customer network

HL7

LIS or middleware



BD Veritor™ Plus Analyzer deployed as a self-contained instance on existing customer infrastructure.

# BD Veritor™ Connectivity Starter Pack

## BD Veritor™ Connectivity Starter Pack + Tests + Analyzer + InfoWifi

### Option 1: Flu + COVID Connectivity Starter Pack

BD Product No.	Quantity	Value Pack includes:	Individual Prod. No.
256139	2	BD Veritor™ System for Rapid Detection of <b>Flu A+B (CLIA-waived)</b> (30 tests per kit)	256045
	2	BD Veritor™ System for <b>SARS-CoV-2 (CLIA-waived)</b> (30 tests per kit)	256092
	1	BD Veritor™ <b>Plus Analyzer</b>	256066
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010

### Option 2: Multi-Assay Connectivity Starter Pack

BD Product No.	Quantity	Value Pack includes:	Individual Prod. No.
256140	2	BD Veritor™ System for Rapid Detection of <b>SARS-CoV-2 &amp; Flu A+B**</b> (30 tests per kit)	256088
	1	BD Veritor™ System for Rapid Detection of <b>RSV (CLIA-waived)</b> (30 tests per kit)	256038
	1	BD Veritor™ System for Rapid Detection of <b>Group A Strep (CLIA-waived)</b> (30 tests per kit)	256040
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010
	1	BD Veritor™ <b>Plus Analyzer</b>	256066

### Option 3: High Volume Triplex Connectivity Pack (SKU 256141)

BD Product No.	Quantity	Value Pack includes:	Individual Prod. No.
256141	3	BD Veritor™ System for Rapid Detection of <b>SARS-CoV-2 &amp; Flu A+B**</b> (30 tests per kit)	256088
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010
	1	BD Veritor™ <b>Plus Analyzer</b>	256066

## Connectivity Expansion Bundle Tests + InfoWifi Module

### Option 1: Triplex Connectivity Bundle

BD Product No.	Quantity	Value pack includes:	BD Product No.
256142	2	BD Veritor™ System for Rapid Detection of <b>SARS-CoV-2 &amp; Flu A+B**</b> (30 tests per kit)	256088
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010

### Option 2: COVID Connectivity Bundle

BD Product No.	Quantity	Value pack includes:	BD Product No.
256143	3	BD Veritor™ System for <b>SARS-CoV-2 (CLIA-waived)</b> (30 tests per kit)	256092
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010

### Option 3: Flu Connectivity Bundle

BD Product No.	Quantity	Value pack includes:	BD Product No.
256144	3	BD Veritor™ System for Rapid Detection of <b>Flu A+B (CLIA-waived)</b> (30 tests per kit)	256045
	1	BD Veritor™ <b>InfoWiFi Module</b>	445010

Order the items above through your distribution partner. BD will drop ship items together directly to you.

# Invest in simplicity

with one system, flexible test menus,  
and streamlined connectivity—without  
added complexity.



Get connected.



Contact your authorized BD distributor or call  
**1.844.823.5433** (say “Veritor” at the prompt)

\*Result processing times for each BD Veritor™ assay are 15 minutes for the SARS-CoV-2 and SARS-CoV-2 & Flu A+B assays, 10 minutes for Flu A+B and RSV assays, and 5 minutes for the Group A Strep assay.

\*\*In the USA, the BD Veritor™ System for Rapid Detection of SARS-CoV-2 & Flu A+B has not been FDA cleared or approved but has been authorized by the FDA under an Emergency Use Authorization for use by authorized laboratories; use by laboratories certified under the CLIA, 42 U.S.C. §263a, that meet requirements to perform moderate, high, or waived complexity tests. The product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

This product has been authorized only for the detection of proteins from SARS-CoV-2, influenza A and influenza B, not for any other viruses or pathogens; and, in the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID 19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. §360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Connectivity agreement will need to be executed prior to EMR integration.

BD Diagnostics, 7 Loveton Circle, Sparks, MD 21152, U.S.  
1.800.638.8663

[bd.com](https://www.bd.com)



BD is a trademark of Becton, Dickinson and Company. Waters is a trademark of Waters Corporation or its affiliates. All other trademarks are the property of their respective owners. © 2026 Waters Corporation. (BD-177429 4848-US-0526 June 2026)