

There's a life behind every test

CLIA-waived testing for Group A Strep

Achieve accurate Group A Strep results in 5 minutes with the portable, easy-to-use BD Veritor™ Plus System.



*EUA Authorized by FDA



Efficient Group A Strep testing

GROUP A STREP

The **BD Veritor™ Plus System** is designed to offer simple Group A Strep testing, with results delivered in just 5 minutes.



Simplifies the testing process

- Easy operation and single-button functionality may help reduce potential manual errors
- Provides intuitive sample processing with pre-filled, unitized tubes



Delivers workflow efficiency

- Adapts easily to your workflow by offering 2 operational modes
 - » **Walk Away:** the test device is inserted into the Analyzer immediately after preparation, allowing staff to multitask while the sample incubates
 - » **Analyze Now:** the test device is inserted after the incubation time is complete, allowing batches of samples to be tested at once



Achieves reliable, rapid results

- Displays easy-to-read digital results for Group A Strep in 5 minutes¹
- Provides clear, objective results by correcting for nonspecific binding and detecting positives not recognized by the unaided eye¹



Provides result traceability

- Downloadable lot number, patient/specimen ID, operator ID entry capability and test records with the BD Veritor™ InfoWiFi module or via the companion ImageMover application
- Offers result printing capabilities via USB port

Reference: 1. BD Veritor™ Plus System for Rapid Detection of Group A Strep, laboratory kit [package insert], 8087675. Franklin Lakes, NJ: Becton, Dickinson and Company

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999 USA
800-638-8663

BDVeritor.com

BD, the BD Logo and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. © 2023 BD. All rights reserved. (BD-43403 2040-US-0323 March 2023)

Offers proven performance vs culture¹

	Sensitivity	Specificity	Results in
Throat swab sample	96.6% (95% CI: 92.4%, 98.6%)	95.5% (95% CI: 93.6%, 96.9%)	5

Enable your staff with the BD Veritor™ Plus eLearning platform

- Assign customized training plans by role, product or learning type
- Simplify learning with 5-minute modules
- Easily track training compliance and competency assessments

Ordering information	Cat no.	Qty
BD Veritor™ Plus System Group A Strep CLIA-waived kit	256040	30 tests
BD Veritor™ Plus Analyzer	256066	1 unit

*EUA Information for the BD Veritor™ SARS-CoV-2 and SARS-CoV-2 & Flu A+B assays: These products have not been FDA cleared or approved; but have been authorized by FDA under EUA for use by authorized laboratories. The BD Veritor™ System for Rapid Detection of SARS-CoV-2 has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; the BD Veritor™ System for Rapid Detection of SARS-CoV-2 & Flu A+B has been authorized only for the detection of proteins from SARS-CoV-2, influenza A and influenza B, not for any other viruses or pathogens; and, the emergency use of these products is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

