


A circular inset image showing a middle-aged couple smiling and looking at each other in a park setting. The man is wearing a light blue checkered shirt and glasses, and the woman is wearing a blue and white striped shirt with lace sleeves.

# Invest in simplicity


Rely on the BD Veritor™ Plus System for portable and rapid point-of-care respiratory testing with results in 15 minutes or less.\*



# Simple means...


## No guessing

Visually-read assays are interpreted subjectively, which can potentially reduce test sensitivity and lead to false results.<sup>1,2</sup>




## No long waits

Delayed treatment decisions negatively impact patient experience.<sup>3</sup>



## No running around

Fixed testing locations away from the patient make testing more complicated and time-consuming.<sup>4</sup>



## No added costs

Committed contracts and monthly maintenance impacts your profitability.<sup>5</sup>



# BD Veritor™ Plus System

Your simple solution



## Digital clarity:

- Unambiguous positive/negative display which may help to reduce result interpretation errors and streamline diagnosis.




## Rapid results:

- Delivers digitally read results in 15 minutes or less,\* supporting timely clinical decisions.



## Easy to use:

- No monthly calibration required, operated with single-button functionality.



## Flexible workflow options:

- Walk Away mode: Built in timer allows you to perform other tasks throughout test development.
- Analyze Now mode: Enables immediate processing to support batch testing.
- Compact, lightweight and portable—able to function without cords.



## Cost-effective:

- No service contracts, maintenance or calibration needed for device lifetime.

Confident, quality care

## Prioritize efficiency and reliability with BD Veritor

The BD Veritor™ Plus System is a point-of-care testing solution for SARS-CoV-2, SARS-CoV-2 & Flu A+B,<sup>†</sup> Flu A+B, Group A Strep, and RSV. It offers:

- The BD Veritor™ Plus Analyzer, a simple, easy-to-use device that delivers unambiguous, rapid and reliable digital results
- A range of assays to aid clinicians in the diagnosis and management of multiple respiratory infections for symptomatic patients so they can take their next steps with more confidence



Tests	SCoV-2 FLU A+B	SARS-CoV-2	FLU	RSV	STREP A
Results in	15 mins <sup>4</sup>	15 mins <sup>5</sup>	10 mins <sup>6</sup>	10 mins <sup>7</sup>	5 mins <sup>8</sup>

Reliable results in a timely manner allow you to quickly diagnose and improve care for your patients<sup>2,3</sup>

RSV=respiratory syncytial virus; SARS-CoV-2=severe acute respiratory syndrome coronavirus  
<sup>†</sup>Refer to the back cover.



## Testing workflow

# Combine easy sample collection and testing workflows for reliable answers you can use

The BD Veritor™ Plus System offers operational flexibility, with 2 modes designed to meet your workflow needs.



### Analyze Now Mode

Batch\*\* multiple tests simultaneously  
(default mode)

### Walk Away Mode

Plug in and the device will time the test  
(2 clicks away)


The BD Veritor™ Plus System also provides simple and efficient workflows that can be incorporated across a variety of clinical settings.

#### Step 1



Collect sample<sup>†</sup>

#### Step 2



Mix swab with  
vial reagent

#### Step 3




Prepare assay  
with sample

#### Step 4



Run test

#### Step 5



Insert cartridge and  
read the results

\*\*Batch testing refers to Triplex and SARS-CoV-2 standalone tests only.

<sup>†</sup>For Group A Strep, add 3 drops of Reagent 1 to the Reagent 2 vial and inoculate the sample for 1 to 2 minutes prior to assay preparation.

# Test. Connect. Report. Print.

The BD Veritor™ Plus System offers several options for connectivity to help you and your team manage and report test results.



## BD Veritor™ InfoWiFi Module

A barcode scanner that captures and links additional information with each result.

## Workflow Services by ImageMover

Connect point-of-care testing to the EHR with streamlined BD Veritor™ results integration.

## BD Informatics

Printing capabilities.

Capability	Analyzer	Analyzer <sup>1</sup> + InfoWiFi <sup>1</sup>	Analyzer <sup>1</sup> + Workflow Services <sup>1</sup>	Analyzer <sup>1</sup> + InfoWiFi <sup>1</sup> + Workflow Services <sup>1</sup>
Rapid, digital results	✓	✓	✓	✓
Identifiable results <sup>5</sup>		✓	✓	✓
Downloadable results		✓	✓	✓
Public health reporting			✓	✓
Automatic real-time results transmission				✓
Integration with EMR/LIS				✓

**Together we can customize an effective and connected reporting solution for your practice.**

EMR=electronic medical records; LIS=laboratory information systems.

<sup>1</sup>Requires a barcoded workflow.


<sup>5</sup>Test result is associated with Patient ID, Operator ID, and Kit Lot information.

# How to order the BD Veritor™ Plus System

Contact your authorized distributor partner or BD to learn more.

For expedited service, say **Veritor** at the prompt. **1.844.823.5433**

BD Product No.	Product Description	Quantity
256066	BD Veritor™ Plus Analyzer	1 unit
256088	BD Veritor™ System for Rapid Detection of SARS-CoV-2 & Flu A+B <sup>+</sup> (EUA)	30 tests
256092	BD Veritor™ System for SARS-CoV-2	30 tests
256045	BD Veritor™ System for Rapid Detection of Flu A+B (CLIA-waived)	30 tests
256041	BD Veritor™ System for Rapid Detection of Flu A+B (moderate complexity)	30 tests
256038	BD Veritor™ System for Rapid Detection of RSV (CLIA-waived)	30 tests
256042	BD Veritor™ System for Rapid Detection of RSV (moderate complexity)	30 tests
256040	BD Veritor™ System for Rapid Detection of Group A Strep (CLIA-waived)	30 tests
445010	BD Veritor™ InfoWiFi Module	1 unit
443907	USB Printer Cable	1 unit



Scan here to view the clinical performance of the BD Veritor™ Plus System.

## References

1. Centers for Disease Control and Prevention. Treatments your healthcare provider might recommend if you are sick. Updated January 13, 2022. Accessed March 18, 2022. <https://www.cdc.gov/coronavirus/2019-ncov/your-health/treatments-for-severe-illness.html>
2. Huang G, Gong T, Wang G, et al. Timely diagnosis and treatment shortens the time to resolution of coronavirus disease (COVID-19) pneumonia and lowers the highest and last CT scores from sequential chest CT. *AJR Am J Roentgenol.* 2020;215(2):367-373.
3. Koonin LM, Patel A. Timely antiviral administration during an influenza pandemic: key components. *Am J Public Health.* 2018;108(Suppl 3):S215-S220.
4. BD Veritor™ System for Rapid Detection of SARS-CoV-2 & Flu A+B. Package insert, 500051910. Becton, Dickinson and Company.
5. BD Veritor™ System for SARS-CoV-2. Package insert. 500075511(01) Becton, Dickinson and Company.
6. BD Veritor™ System for Rapid Detection of Flu A+B. Package insert, 8087667. Becton, Dickinson and Company.
7. BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV). Package insert, 8086098. Becton, Dickinson and Company.
8. BD Veritor™ System for Rapid Detection of Group A Strep. Package insert, 8087675. Becton, Dickinson and Company.

The BD Veritor™ Plus System

## Training solutions tailored to support your team and instill confidence

**BD offers support and resources based on your needs, so you and your team can feel confident using the BD Veritor™ Plus System.**

Whether you are new to BD Veritor or a seasoned user, BD provides a variety of learning options to help you and your team confidently administer trusted, reliable point-of-care testing.



Scan to **access training modules** on the BD Veritor™ website



Contact us to request a **demo**



Learn more at **bdveritor.com**

\*Result processing times for each BD Veritor™ assay are 15 minutes for the SARS-CoV-2 and SARS-CoV-2 & Flu A+B assays, 10 minutes for Flu A+B and RSV assays, and 5 minutes for the Group A Strep assay.

\*In the USA, the BD Veritor™ System for Rapid Detection of SARS-CoV-2 & Flu A+B has not been FDA cleared or approved but has been authorized by the FDA under an Emergency Use Authorization for use by authorized laboratories; use by laboratories certified under the CLIA, 42 U.S.C. §263a, that meet requirements to perform moderate, high, or waived complexity tests. The product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

This product has been authorized only for the detection of proteins from SARS-CoV-2, influenza A and influenza B, not for any other viruses or pathogens; and, in the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID 19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. §360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152, U.S.  
1.800.638.8663

**bd.com**



BD, the BD Logo, BD Synapsys and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved. (BD-158014 4485-US-1125 November 2025)