


# COVID-19 TESTING IS HERE

Keeping our school safe is job #1 and that's why we've partnered with BD, an industry leader in rapid antigen testing, to provide onsite COVID-19 diagnostic testing right here in our school:



**Comfortable  
nasal swab**



**Results in  
15 minutes**



**Proven, digitally  
read results**



If you have any questions about COVID-19 testing in our school, please contact:

Should you have questions about your COVID-19 testing program, please contact **1-888-4BDCOVID** and **veritorcares@bd.com**

\* Information about the BD Veritor™ System for Rapid Detection of SARS-CoV-2:

- This product has not been FDA cleared or approved; but has been authorized by FDA under an EUA for use by authorized laboratories
- This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and
- This product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 is a chromatographic digital immunoassay intended for the direct and qualitative detection of SARS-CoV-2 nucleocapsid antigens in nasal swabs from individuals who are suspected of COVID-19 by their healthcare provider within the first five days of the onset of symptoms. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet the requirements to perform moderate, high or waived complexity tests. This test is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

**BDVeritor.com**

BD, the BD Logo, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates.  
All other trademarks are the property of their respective owners. © 2021 BD. All rights reserved. BD-868-US-0321 (03/21)

