

There's a life behind every test

Point-of-care testing for Group A Strep

Achieve accurate Group A Strep results in 5 minutes with the portable, easy-to-use BD Veritor™ Plus System.



Efficient Group A Strep testing

GROUP A STREP

The **BD Veritor™ Plus System** is designed to offer simple Group A Strep testing, with results delivered in just 5 minutes.



Simplifies the testing process

- Easy operation and single-button functionality may help reduce potential manual errors
- Provides intuitive sample processing with prefilled, unitized tubes



Delivers workflow efficiency

- Adapts easily to your workflow by offering 2 operational modes
 - » **Walk Away:** the test device is inserted into the Analyzer immediately after preparation, allowing staff to multitask while the sample incubates
 - » **Analyze Now:** the test device is inserted after the incubation time is complete, allowing batches of samples to be tested at once



Achieves reliable, rapid results

- Displays easy-to-read digital results for Group A Strep in 5 minutes¹
- Provides clear, objective results by correcting for nonspecific binding and detecting positives not recognized by the unaided eye¹




Provides result traceability

- Downloadable lot number, patient/specimen ID, operator ID entry capability and test records with the BD Veritor™ InfoWiFi Module or via the companion ImageMover application
- Offers result printing capabilities via USB port

Reference: 1. BD Veritor™ System for Rapid Detection of Group A Strep, laboratory kit [package insert], 8087675. Franklin Lakes, NJ: Becton, Dickinson and Company

BD Life Sciences, 2100 Derry Road West, Mississauga, ON L5N 0B3 CA
866-979-9408

Offers proven performance vs culture¹

	Sensitivity	Specificity	Results in
Throat swab sample	96.6% (95% CI: 92.4%, 98.6%)	95.5% (95% CI: 93.6%, 96.9%)	

Enable your staff with the BD Veritor™ Plus System eLearning platform

- Assign customized training plans by role, product or learning type
- Simplify learning with 5-minute modules
- Easily track training compliance and competency assessments

Ordering information	Cat no.	Qty
BD Veritor™ Plus System Group A Strep kit	256040	30 tests
BD Veritor™ Plus Analyzer	256066	1 unit

¹BD Veritor™ System for Rapid Detection of SARS-CoV-2:

- This test has been authorized for sale in Canada by Health Canada under Interim Order;
- This test is only authorized for the duration of the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19, unless the authorization is terminated or revoked sooner;
- This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.

BDVeritor.com

BD, the BD Logo and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved.

